


Best 10 Best 10 Best
Best 10 Best 10 Best
Best 10 Best 10 Best


WHICH
10 SECONDS
WOULD YOU
CHOOSE?




 **Mortara**

The

last 10 s




econds?




Mortara's Best 10 algorithm evaluates ECG noise to determine the best quality 10 seconds of data recorded by the ELI™ electrocardiograph. So, no matter when you press the ECG button, you get the best quality ECG for review - even if your patient moved in the last 10 seconds.

or the

BEST 10



seconds?



Fewer Repeat ECGs
Less Time Wasted
Better Quality Recordings

Mortara is committed to simplifying
everyday ECG
acquisition.





MORTARA INSTRUMENT, INC.

7865 North 86th Street

Milwaukee, WI 53224

U.S.A.

Tel: 414.354.1600

Tel: 800.231.7437

Service: 888.MORTARA

Fax: 414.354.4760

ELI™ is a trademark of Mortara Instrument, Inc.

