



see the
TRUTH
in ECG waveforms



 **Mortara**

Diagnostic Quality ECG


Is your electroca



From H. Burri et al. "Simulation of anteroseptal myocardial infarction by electrocardiographic filters", *Journal of Electrocardiology* 39 (2006) 253-258

radiograph

lying to you?




If its filter setting is above .05 Hz,

it probably is!

Non-diagnostic settings produce nice, clean tracings, but may also distort waveform morphology, potentially leading to false-negative and false-positive diagnoses.


Diagnostic quality ECG matters!

Don't be lied to. Choose a diagnostic filter setting and see the truth in ECG waveforms.



radiograph

lying to you?





MORTARA INSTRUMENT, INC.

7865 North 86th Street

Milwaukee, WI 53224

U.S.A.

Tel: 414.354.1600

Tel: 800.231.7437

Service: 888.MORTARA

Fax: 414.354.4760

